

ADVANCED ENDODONTIC SYSTEMS®

ERGO V DRAWER/DOOR SYSTEMS

The New Ergo V

The newest addition to ASI's *Designer Series* line of Advanced Endodontic Systems, the Ergo V, takes customization to a new level. The Ergo V features a unique, angled shape to provide improved visibility of instrument control panels and increased worktop space.

The Ergo V can be ordered left or right facing to provide ergonomic delivery from either side of the patient. By adding the Assistant's Instrument Package (air/water syringe, saliva ejector, high volume evacuator, and solids trap), the system becomes an ideal rear delivery unit.

New Features

The Ergo V incorporates the latest advancements in design technology. New handpiece holders improve precision and safety during activation and provide exceptional access to your instruments. Stainless steel components permit worry free disinfection and enhanced reliability. The new chassis design, featuring fan cooling and a flip-down control panel, facilitates servicing and upgrades to the system. New custom colors and bonded quartz worktops add impact to office decor.




DESIGNER SERIES


*Ergo V Drawer/Door Cart
Model 2250E*

*Shown:
Right Facing,
Rear Delivery System*

*Also available in
left handed configurations.*



◀ *Rear Delivery
Ergo V with
Assistant's Package*



▶ *Side Delivery
Ergo V with
Assistant's Cart*

“ASI worked to create a modern version of our delivery system, integrating a computer system, and integrating all of our instruments with only one foot control...ASI is a progressive dental company!”

*Dr. Gabriel Tagher
Cape Cod, Massachusetts*

ERGO V™ DRAWER/DOOR SYSTEM

IMPROVED DESIGN FEATURES

► *Personalized Integration*

Six instrument panels are available for the ultimate in instrument integration. Ultrasonic units, obturation systems, electric motors, and apex finders are conveniently located and individually configured for fingertip operation. All functions operate from a single foot control.

► *Capability*

Up to three pneumatic 4-hole handpiece connectors are included with an air/water syringe and an Endodontic air only syringe with precision regulated low pressure.

► *Functionality*

A dual clean closed water system is a standard feature located conveniently and unobtrusively behind the hinged door. The four drawer storage system can be specified left or right facing. The system also features a built-in nook for additional instrument placement.

► *Upgradeable*

The non-porous worktop area can be enhanced with optional swivel trays or CPU, monitor, and keyboard mounts. The large caster base provides stability for the optional mounting of microscopes and prevents tubing drag.



*Shown:
Left Facing,
Side Delivery System*

*Also available in
left handed configurations.*

*Sample configurations shown with optional instruments.
Systems will be custom integrated with the brands/models
of instruments you specify and in the order of placement
you request. All systems can be equipped with optional
computer monitor, keyboard and CPU mounts.*

Prices and specifications subject to change without notice and do not include shipping and handling.

